

The logo for Precision Instruments Corp features the company name in a white, sans-serif font. The text is split across two overlapping horizontal bars: an orange bar on the left and a black bar on the right. The background of the entire page is white, with a large, vibrant blue curved shape on the left side that partially overlaps the orange bar.


Precision Instruments Corp

Surgical Instruments


Atrauma-Gefaßklemmen

Atrauma Artery Forcep and clamps


Pinzas Clamps Atrauma




Toothing according to Cooley:
 Very narrow jaws for application in extremely limited operational field. Particularly approved for cardiovascular surgery and pediatry. Firm grip without squeezing.




Cooley
 1-464-16 = 16 cm / 6 1/4"
 1-464-20 = 20 cm / 8"



Cooley
 1-468-17
 17 cm / 6 3/4"



Cooley
 1-469-20
 20 cm / 8"



Cooley
 1-471-16
 16 cm / 6 1/4"

Atrauma-Gefaßklemmen

Atrauma Artery Forcep and clamps


Pinzas Clamps Atrauma



Cooley
1-472-15
15 cm / 6"



Cooley
1-472-16
15,5 cm / 6 1/8"



Cooley
1-472-17
17 cm / 6 3/4"



Cooley
1-473-17
17 cm / 6 3/4"



Cooley
1-474-17
17 cm / 6 3/4"



Cooley
1-476-17
17 cm / 6 3/4"

Atrauma-Gefaßklemmen

Atrauma Artery Forcep and clamps

Pinzas Clamps Atrauma



Cooley-Derra
1-478-16 - 1-478-18



1-478-16
16 cm / 6 1/4"



1-478-17
17 cm / 6 3/4"



1-478-18
18 cm / 7"



Cooley
1-477-16 = 16 cm / 6 1/4"
1-477-17 = 17 cm / 6 3/4"



1-477-16



1-477-17



Beck
1-479-17
17 cm / 6 3/4"



1/2

