


The logo for Precision Instruments Corp features the company name in a white, sans-serif font. The text is split across two horizontal bars: an orange bar on the left and a black bar on the right. The background of the entire page is white, with a large blue curved shape on the left side.


Precision Instruments Corp

Surgical Instruments


Plazentar - und Abortuszangen Placenta And Ovum Forceps Pinzas para aborto y placenta



Winter
28 cm/11"




straight	Fig.	curved
1-344-01	1	1-345-01
1-344-02	2	1-345-02



Kranz
1-360-29
29 cm/11 1/2"

1-361-29



The holding part has been shifted to the inside of the blades. The uterine mucosa will therefore not be gripped by the edge of the jaws.



Corey (Greenhalg)
1-363-24
24 cm/9 1/2"



Cuzzi
30 cm/12"



1/4
toothed back



14 1/2 mm
1-382-01



18 mm
1-382-02



21 mm
1-382-03

blunt