

The logo for Precision Instruments Corp features the company name in a white, sans-serif font. The text is split across two overlapping horizontal bars: an orange bar on the left and a black bar on the right. The background of the entire page is white, with a large, vibrant blue curved shape on the left side that partially overlaps the orange bar.

Precision Instruments Corp

Surgical Instruments

INDEX

1. Laryngoskope
Laryngoscopes 3 to 3 page
Laringoscopios


2. Katheterzangen
Catheter Forceps 4 to 4 page
Pinza para cateteres

Laryngoskope Laryngoscopes Laringoscopios



McIntosh
1-380-00

19-380-00
Laryngoscope Set,
consisting of:
1-392-00 Battery Handle (without batteries)
1-383-01 McIntosh Blade, Fig. 1
1-383-02 ditto, Fig. 2
1-383-03 ditto, Fig. 3
1-398-03 Case for 3 curved blades



1-383-00

1-383-01



1-383-02



1-383-03




1-383-04

McIntosh Laryngoscope Blades
1-383-00 Fig. 0, working length 55 mm, baby's
1-383-01 Fig. 1, working length 75 mm, children's
1-383-02 Fig. 2, working length 90 mm, adolescents'
1-383-03 Fig. 3, working length 110 mm, women's
1-383-04 Fig. 4, working length 135 mm, men's

Katheterzangen

Catheter Forceps

Pinza para cateteres



1-250-02

1-301-20 – 1-303-25

Oxygen Insufflating Spectacles with
2 tubing tips
Catheter Introducing Forceps